

# ACCESSING MAMMALIAN DRUG METABOLITES USING POLYCYPS® ENZYMES


Sandie Lai, William Hodds, Vincent Poon, Kinga Nytko, Aksana Khan, Emily Hopkins, Antonio de Riso, Richard Phipps, Jonathan Steele, Stephen Wrigley, Liam Evans, Frank Scheffler, Julia Shanu-Wilson

Hypha Discovery, 957 Buckingham Avenue, Slough, SL1 4NL, UK. Email: julia.shanuwilson@hyphadiscovery.co.uk



**Abstract:** This poster illustrates the application of a new biocatalysis kit, PolyCYPs®, to enable scalable synthesis of CYP-derived metabolites of drugs. The PolyCYPs platform is comprised of a set of recombinant cytochrome P450 enzymes and redox partners cloned from some of the talented actinomycetes in Hypha's biotransformation panel. Enzymes in the kit catalyze the oxidation of a wide variety of substrate types to generate multiple mammalian and microbial-derived CYP metabolites. The poster features application of selected PolyCYPs isoforms to produce oxidized human metabolites of drugs. Further, the utility of PolyCYPs enzymes for introducing oxygen into a drug candidate as part of a late stage functionalization program is illustrated, in which derivatives can be generated in parallel which may possess superior properties such as improved metabolic stability and LLE (lipophilic ligand efficiency), and for exploration of structure-activity relationships.

PolyCYPs enzymes are proprietary recombinant microbial CYPs expressed in *E. coli* together with ferredoxin / ferredoxin reductase redox partners. They are derived from a subset of Hypha's actinomycetes talented at oxidizing drugs.



## Lyophilized and vacuum-sealed reagents included in PolyCYPs screening kits:

PolyCYPs enzymes (**blue lid**) - PolyCYPs include 6, 152, 166, 168, 194, 196, 217, 235

Co-factor vial (**green lid**) - NADP+, G6P, G6PDH

Formulant (**red lid**) - for highly lipophilic substrates

Control compound (**silver lid**) - bosentan

## Process summary for PolyCYPs screen

Lyophilised PolyCYPs enzymes reconstituted in water according to PolyCYPs kit instructions



Drugs dosed in DMSO at 0.1 mg/ml in 0.4% final DMSO concentration (up to 2% possible)

Lyophilised co-factors reconstituted in water and added to the PolyCYPs and drug substrate to initiate the reaction

Shaken overnight (16-20 hrs) in 24 square well plate at 27°C

Reaction stopped with equal volume of acetonitrile, centrifuged and supernatant analysed by LC-MS

Scan QR code for a pdf of the poster




## Generation of human CYP-derived metabolites

### Results

#### Step 1: Multiple drugs (20) screened using 8 PolyCYPs enzymes

> 10% conversion of 18 drugs to at least one metabolite



#### Meloxicam

Human metabolism: CYP2C9 (CYP3A4)

5'-hydroxymethyl meloxicam (M5)

Major human *in vivo* metabolite

Max 66% conversion

5'-carboxy meloxicam (M5)


Major human *in vivo* secondary metabolite (low turnover in HLMs)

7.4% conversion

#### Ritonavir

Human metabolism: CYP3A4 & CYP2D6

7 enzymes giving 28-90% conversion



Human metabolite

$\epsilon$ -HO-N-methylisopropylthiazolylmethyl-ritonavir (M2)


#### Step 2: Scale-up of reaction of interest

Multiple scale-up vials, or fresh enzyme prep, can be used to generate multi mg amounts of metabolites. For larger quantities, scale-up is achieved using a recombinant *Streptomyces* clone expressing the PolyCYP undertaking the biotransformation of interest, or the originating wild type strain.

#### Ritonavir scale-up

Hydroxylated derivatives of ritonavir produced in an unoptimized screening reaction with PolyCYP 168

Hydroxylated derivatives of ritonavir produced in a bioreactor by a *Streptomyces* clone containing PolyCYP 168



## Conclusions

- PolyCYPs enzymes can be used to generate mg quantities of human metabolites and other hydroxylated derivatives of drugs for MetID & biological evaluation.
- Recombinant streptomycete clones expressing PolyCYP enzymes enable large scale production of synthetically challenging metabolites / derivatives.

## Late stage oxidation of drug compounds

- Previously we generated several derivatives of one of AstraZeneca's drug leads by microbial whole cell biotransformation.
- Analogues were characterised/quantified by 2D NMR/qNMR, and found to be derived from oxidation of the cyclohexane moiety, together with desmethyl and benzylic hydroxylated derivatives.
- Subsequently we investigated the ability of 5 PolyCYPs enzymes (6, 14, 152, 166, 168) to generate the same oxidized derivatives.

## Results

- All five PolyCYPs enzymes used generated some hydroxylated derivatives of the parent drug (mock structures shown due to confidentiality) although PolyCYPs 152 and 166 generated the best array for this parent compound.
- PolyCYP 166 gave a 30% conversion to 4 monohydroxylated derivatives.
- PolyCYP 152 gave a 58% conversion to produce 5 derivatives that matched with all monohydroxylated compounds previously observed.

